MTA REPAIR HP

PORTUGUÊS

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ENGLISH

Bioceramic reparative cement with high plasticity

INTRODUCTION

MTA (Mineral Trioxide Aggregate) has its established use in Endodontics due to the clinical results proven through numerous scientific studies. However, due to the physical properties inherent to the powder of the product, the "sandy" consistency hinders manipulation and transport of the material to the site of repair. To improve these characteristics, ANGELUS® has developed a new formulation called MTA REPAIR HP – "High Plasticity" MTA. This new formula maintains all the chemical and biological properties of the original MTA, which guarantees the success of the treatment, but changes its physical properties of manipulation. The result is a product with greater plasticity, facilitating manipulation and insertion in the dental cavity.

DEFINITION

Restorative cement ready for use.

MTA REPAIR HP is an endodontic restorative cement with high plasticity, composed of mineral oxides in the form of fine hydrophilic particles. It is indicated for cases of root perforation (canal and furcation), iatrogenic or by caries, root perforation by internal resorption, retrofilling, direct pulp capping, pulpotomy, apexigenesis, and apexification.

Advantages over other materials:

- New formula whose particle size after hydration allows for easy manipulation and insertion into the dental

- cavity;
- Addition of a radiopacifier (CaWO_4) which does not cause staining of the root or dental crown;
 - Low solubility;
 - Excellent marginal sealing capacity through controlled setting expansion which impedes the migration of microorganisms and fluids to the inside of the root canal;
 - Excellent biological sealing of root perforation (channel and furcation) to induce formation of periradicular cement;
 - Promotes the formation of a dentin bridge when used in pulp capping;
 - Allows for use in wet medium without change of its properties.

COMPOSITION

- Powder: Tricalcium silicate $3\text{CaO}\cdot\text{SiO}_2$; Dicalcium silicate $2\text{CaO}\cdot\text{SiO}_2$; Tricalcium aluminate $3\text{CaO}\cdot\text{Al}_2\text{O}_3$; Calcium oxide CaO ; Calcium Tungstate CaWO_4 ;
- Liquid: Water and Plasticizer.

INDICATIONS

1. Treatment of perforations of root canal and furcation caused iatrogenically or by caries lesion;
2. Via canal treatment of root perforation due to internal resorption;
3. Surgical treatment of root perforation due to internal resorption;
4. Periapical surgery with reverse filling;
5. Pulp capping;
6. Pulpotomy (removal of the affected crown portion of the pulp preserving the vitality and function of the remaining radicular portion);
7. Apexogenesis (induction of root development in vital teeth with an inflamed coronal pulp);
8. Apexification (induction of formation of a mineralized barrier at the root tip of young permanent teeth with incomplete root development and a necrotic pulp).

PROPERTIES

- Setting time: MTA REPAIR HP solidifies when kept in a wet environment after spatulation. Initial setting time is approximately 15 minutes;
- Radiopacity: Nearly matches that of gutta-percha. More radiopaque than dentine and bone;
- Absence of dental discoloration due to the CaWO_4 radiopacifier used.

TECHNIQUE OF USE

1. Sterilize the materials to be used in spatulation, insertion, and condensation of the MTA REPAIR HP;

2. Dispense the content of 1 package of MTA REPAIR HP and 2 drops of the liquid on the glass tray;
3. Spatulate for 40 sec. until the powder and liquid are completely homogenized. The cement obtained is similar to a modeling compound;
4. Take the MTA REPAIR HP to the desired location with an ANGELUS® MTA APPLICATOR, condensing it with suitable instruments for this purpose.

PRECAUTIONS

- Use eye protection, mask and gloves when handling MTA REPAIR HP. In case of eye or skin contact, wash abundantly with water;
- Open the package of MTA REPAIR HP only immediately before use. It is very sensitive to moisture;
- Do not use MTA REPAIR HP to fill a root canal. Its viscosity is inadequate for this procedure;
- Do not use MTA REPAIR HP on the tooth surface that is in contact with the gingival sulcus, because disintegration of the cement will occur;
- Only use MTA REPAIR HP after remission of acute signs and symptoms of the endodontic disease. The acidic pH of endodontically compromised sites (lesions) prevents its setting reaction;
- When applying the MTA REPAIR HP, avoid letting it overflow to regions beyond the dental cavity, because although even if the excess is reabsorbed, it can slow healing.

ATTENTION: The information provided in this manual is based on laboratory and clinical studies. The successful use of MTA depends on a correct diagnosis, the operative technique, the condition of the treated tooth and the general health of the patient. This product must be used according to this manual.

STORAGE

Store the product in a dry place protected from light.

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