# FIBREX-LAB SISTEMA SYSTEM

## PORTUGUÊS

### INTRODUÇÃO

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## ENGLISH

### INTRODUCTION

Glass fiber was developed as a solution to combine strength and lightness. Although limited, its use in Dentistry started in the 1960's. With further advances in silane chemistry and adhesion technology, the 80's marked a true glass fiber revolution in Prosthodontics. The esthetic characteristics of glass fiber along with its most desirable mechanical properties allow immediate, highly conservative, durable and low cost restorative and preventive dental work.

### DESCRIPTION

The FIBREX-LAB SYSTEM is formed by glass fibers impregnated with light-cured composite resin and used for the structural reinforcement of composite resin and ceromer bonded and fixed bridges.

### INDICATIONS

Structural reinforcement of indirect composite resin and ceromer restorations such as inlays, onlays, single crowns, bonded and fixed bridges up to three units. To be used in dental labs only.

### COMPOSITION

- Glass Fibers: Uni, multidirectional and braided fiber glass.
- Impregnated Resin: Bis-GMA resin, urethane dimethacrylate, barium ceramic glass, highly dispersed silicon dioxide, catalysts, pigments.
- Adhesive Resin C: Bis-GMA resin, urethane dimethacrylate, catalysts, pigments.
- Adhesive Resin F: Bis-GMA resin, urethane dimethacrylate, barium ceramic glass, highly dispersed silicon dioxide, catalysts, pigments.
- Bonding agent: silane solution in alcohol.

### PRECAUTIONS

- FIBREX-LAB SYSTEM components must be protected from light, especially after opening the package.
- When exposed to environmental light, working time of all components that present light-cured composite resin is 4 (four) minutes. These components are FIBREX MEDIAL, JUNCIONAL, CORONAL, RIBBON and ADHESIVES RESIN C and F.
- FIBREX-LAB reinforcement fibers have mechanical properties which will withstand masticatory loads when used according to the indications in adequately prepared teeth. It is also necessary to follow the instructions of each manufacturer of composite resin or ceromer used for the restoration.
- The correct diameter of fibers for each indication must be strictly followed according to the directions.
- The fibers must be covered by the restorative composite resin or ceromer in all areas of the restoration.
- Restorations fabricated with the FIBREX-LAB SYSTEM must be cemented with an adhesive technique and with a resin cement.
- Gloves must be used when handling components of FIBREX-LAB SYSTEM. Contact with non-polymerized material and fiber glass powder may cause skin irritation.
- Finishing and polishing of polymerized components of FIBREX-LAB generate fiber glass particles. It is strongly advised to wear eye protection and, to avoid inhaling, a mask and a suction device are necessary.
- Keep FIBREX-LAB kit in the refrigerator. Do not leave it near laboratory furnaces or other sources of excessive heat, which may cure some of its components.
- To be used in aseptic conditions.

### ADVANTAGES

- Esthetics: as a non-metallic material, its translucency provides great esthetic results.
- Uniform distribution of masticatory stresses.
- Lightness: allows fabrication of slim restorations.
- Lower lab fees compared to ceramic restorations.
- Allows repair in the mouth.
- Compatible to all other composite resin/ceromer lab systems.
- Easy handling: adequate fiber dimensions for all situations.
- Practical: glass fibers previously impregnated with light-cured composite resin.
- Biocompatible.

## DIRECTIONS FOR USE - FIXED BRIDGE WITH INLAY/ONLAY ABUTMENTS

To see more instructions for use: [www.angelus.ind.br](http://www.angelus.ind.br)

- With the special stone working cast ready, apply a die spacer. Bridges are made without sectioning the die.
- Manufacture an acrylic resin structure with extensions for the abutment teeth and a pontic. Check dimensions at the connections between abutments and pontic (ideally a minimum of 3mm for height and breadth) and distance between pontic and working cast (1 mm - detail).
- Place the structure on the duplicated model; envelop it completely with high consistency silicone (impression material) leaving the occlusal area exposed.
- Cut the excess of silicone with a scalpel, exposing more of the occlusal, facial and lingual areas of the acrylic resin structure. Make two dents on the silicone - one oval (a) and another round (b). With a carbide bur, flatten the abutment teeth on the duplicated model by cutting out their cusp tips. Paint a layer of release agent on the model, acrylic resin structure (detail) and surrounding silicone.
- Cover the model, acrylic resin structure and silicone with hot glue (dispensed by a glue gun); wait for the glue to set and cool.
- Remove the glue mold and the structure from the duplicated model; with a small round bur, make holes in areas that approximately correspond to cusp tips of the abutment and pontic teeth; with a disc, make grooves on the silicone (to later allow excess of fiber to seep out).
- Apply release agent on the duplicated model.
- Cut two pieces of FIBREX-JUNCIONAL: one, according to the size of the acrylic resin structure, should be placed in the glue mold; the other, between the abutment teeth, on the silicone. Apply ADHESIVE RESIN C on the abutment teeth (model).
- Cut and place a segment of FIBREX-MEDIAL along the extensions of the teeth preparations (arrow a). Place a shorter segment in the pontic area (arrow b).
- Place the glue mold with the fibers back on the duplicated model; now the glue mold - fiber - die - silicone assembly is taken to the fiber pressing and light-curing equipment to be processed in vacuum for 2 minutes. Light-cure for 6 minutes.
- Remove the glue mold leaving the fiber structure on the model and light-cure for 5 more minutes.
- Remove the excess of fiber with double-faced diamond disks or fine-cut tungsten carbide burs (fiber structure with excess - detail).
- Perform finishing with a diamond-impregnated rubber wheel and fine-cut tungsten carbide burs. Check dimensions: connection sites between abutments and pontic should have a minimum of 3mm for height and breadth; inside the preparations, these dimensions should be at least 0.5mm (detail).
- Place the fiber structure on the working model and check the interocclusal space; ideally 1 to 2mm for composite resin application.
- With 50 to 60um aluminum oxide abrasive and pressure of 1 bar (15 pounds), sandblast the fiber structure. Clean it with steam or with ultrasound in an alcohol bath. Dry it off with clean air.
- Paint a layer of bonding agent (silane) on the fiber structure, let it dry for 1 minute and gently apply air free of oil and humidity.
- Paint release agent on the working cast, place the structure on the cast and apply a thin layer of ADHESIVE RESIN F on the surface where the restorative composite resin will be applied.
- Light-cure for 3 minutes.
- Fiber structure ready for application of restorative resin.
- Final occlusal view of the composite resin bridge.

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**ATTENTION:** This product must be used according to the instructions described in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use.

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