3D Moldex/Bite and Moldex

Instructions for use

Choose your language

- | | |
|-----------------------|-----------------------|
| PT Português | SK Slovensky |
| EN English | PL Polski |
| ES Español | ET Eesti |
| FR Français | LT Lietuviškai |
| IT Italiano | LV Latviski |
| RO Română | CS Český |
| NL Nederlands | GA Gaeilge |
| DE Deutsch | HU Magyar |
| SV Svenska | TR Türkçe |
| DA Danish | MT Malti |
| FI Suomi | BG Български |
| HR Hrvatski | RU Русский |
| SL Slovenščina | EL Ελληνικά |

Português

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English

MOLDEX - DISPOSABLE PARTIAL TRAY FOR DOUBLE IMPRESSION AND BITE REGISTRATION

3D MOLDEX - DISPOSABLE ADJUSTABLE TRAY FOR DOUBLE-ARCH IMPRESSION (PARTIAL OR TOTAL)

3D MOLDEX BITE - Tray for bite registration

DEFINITION

Disposable plastic double trays.

INDICATION

For the simultaneous impression of upper and lower arches and bite registration.

COMPOSITION

Frame..... 100% Polyethylene

Mesh 100% Polypropylene

TECHNIQUE OF USE - MOLDEX AND 3D MOLDEX

A. ONE-STEP IMPRESSION

1. Select the tray and adjust in accordance with the size of the patient's arch, as in the illustrations. Ask the patient to bite down and check for the presence of any interference during closure. Ask the patient to bite down several times on the position marked for impression to practice correct closure;

2. Using a pencil, mark the reference of the tray edge on the patient's teeth to ensure insertion in the correct position;

3. Fill with heavy body impression material (high viscosity) on both sides of the tray (ANGELUS® 3D ADDITION SILICONE);

4. Inject the light body impression material (low viscosity) on the critical areas of the preparations to be copied (ANGELUS® 3D ADDITION SILICONE);

5. Introduce the tray into the mouth, check its position against the pencil mark and ask the patient to bite down;

6. Once the impression material has set, ask the patient to release the bite and remove the tray from the mouth;

7. Cast the mold with plaster (see C. Instructions for casting mold).

ATTENTION: For successful double impression, correct indication and training in the technique are required.

B. TWO-STEP IMPRESSION

1. Select the tray and adjust in accordance with the size of the patient's arch, as in the illustrations. Ask the patient to bite down and check for the presence of any interference during closure. Ask the patient to bite down several times on the position marked for impression to practice correct closure;

2. Using a pencil, mark the reference of the tray edge on the patient's teeth to ensure insertion in the correct position;

3. Fill with heavy body material (high viscosity)

- (ANGELUS® 3D ADDITION SILICONE) - both sides of the tray covering it with a plastic film. This will create relief areas for the second stage of the impression;

4. Introduce the tray into the mouth, check its position against the pencil mark and ask the patient to bite down;

5. Once the impression material has set, ask the patient to release the bite and remove the tray from the mouth;

6. Remove the plastic film from the impression and inject the light body material (low viscosity) into the area previously in relief (ANGELUS® 3D ADDITION SILICONE);

7. Introduce the tray into the mouth again and ask the patient to bite down;

8. Once the light body material has set, ask the patient to release the bite and remove the tray from the mouth.

C. INSTRUCTIONS FOR POURING THE IMPRESSIONS

1. Remove the excess molding material from the edge of the tray;

2. Cast the side of the mold with the dental preparations using plaster. Wait until set, remove the model and die cut it. This first cast will be used as the main working model for manufacturing the prosthetic part (waxing and sealing of the cervical edges of the preparation);

3. Cast the side of the mold with the dental preparations again. Wait for the plaster to set. Check whether there was any overflow of plaster onto the opposite arch. If there was, remove it and apply liquid Vaseline to these points;

4. Keep the second model in the mold with the dental preparations and cast the opposite side;

IMPORTANT: Keep the models in position, without removing them from the impression.

5. With the two models in the mold, fix them with plaster to a simple articulator;

6. Separate the templates and remove the mold (tray and impression material). These templates assembled on the articulator will be used for checking the proximal and occlusal contacts of the prosthesis.

TECHNIQUE OF USE - 3D MOLDEX BITE

1. Test the tray to check for interference;

2. Insert the bite registration material selected on both sides of the tray;

3. Bring it to the patient's mouth, asking them to bite down, holding it until the material sets;

4. Remove it from the mouth with a single movement after the patient releases their bite.

WARNINGS

- Check the mold carefully. The impression should be repeated if the support of the tray interferes in the mold;

- Only use rigid impression material (adding-cure or condensation silicons and polyether);

- Patients with little space for inserting the tray,

due to the presence of third molars, tongue or a bulky torus, may be contraindicated for the double impression technique.

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