Instructions for use

0514102016

Choose your language

- | | |
|-----------------------|-----------------------|
| PT Português | SK Slovensky |
| EN English | PL Polski |
| ES Español | ET Eesti |
| FR Français | LT Lietuviškai |
| IT Italiano | LV Latviski |
| RO Română | CS Český |
| NL Nederlands | GA Gaeilge |
| DE Deutsch | HU Magyar |
| SV Svenska | TR Türkçe |
| DA Danish | MT Malti |
| FI Suomi | BG Български |
| HR Hrvatski | RU Русский |
| SL Slovenščina | EL Ελληνικά |

Português

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Régua calibradora de cones de guta-percha

INTRODUÇÃO

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English

Verify instructions for use.

Callibrating ruler for Gutta-percha cones

INTRODUCTION

This cannot be stated in regard to gutta-percha cones. Due to their plasticity, they may present alterations. Besides the calibration of standardized cones, INTERMEDIUM prepares master cones of intermediate size. This is a great advantage when a standardized cone does not adapt to the apical diameter of root canals.

INDICATIONS

Calibration of standardized gutta-percha cones.

Preparation of intermediate sizes of cones, with in between standard sizes. For example, a size 27.5 cone - between 25 and 30.

COMPOSITION

High performance polymer and metal inserts.

PROPERTIES

Perfect calibration (checking and corrections) of standardized gutta-percha cones, according to the last file used.

Preparation of intermediate sizes of gutta-percha cones, not available in the market.

TECHNIQUE FOR USE

A. Cone calibration

1. Select the master cone according to size of last file used for root canal instrumentation;
2. Place the cone with slight pressure in the perforation that corresponds to size of last file (numbers are written in the back side of the ruler);
3. Verify three possible situations:

The cone enters the perforation and tip reaches evenly the front side surface of ruler: this cone does not need calibration.

The cone enters the perforation but tip does not reach the end of the perforation at the front side surface of ruler: this cone is too large; select another.

The cone passes through the perforation (fig. 1). Cut the tip excess with a new blade (fig. 2).

IMPORTANT: After calibration, margin for error is 0.02 mm.

B. Preparation of intermediate sizes of cones

1. Select the master cone and calibrate it according to the previous instructions;

- Place this cone in the horizontal groove at the lower part of the front side of ruler (fig. 3);
- With a new blade in the vertical groove of the ruler, cut the cone tip. Due to cone taper, the new tip has a larger diameter now. With this procedure, a size 25 cone is converted into a size 27.5 (fig. 4). Cone adequately cut (fig. 5).

IMPORTANT: For sizes 60, 70 and 80, two continuous cuts are necessary to obtain the intermediate sizes.

PRECAUTIONS

Before sterilization, wash the ruler with a neutral detergent and a non-abrasive sponge.

Sterilize it in autoclave (121°C and pressure of 1 atm for 30 minutes) or dry heat (170°C for 2 hours).

Use only new blades to cut the cones.

Chemicals for radiograph processing may stain **INTERMEDIUM**.

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