Instructions for use

0605102016

Choose your language

- | | |
|-----------------------|-----------------------|
| PT Português | SK Slovensky |
| EN English | PL Polski |
| ES Español | ET Eesti |
| FR Français | LT Lietuviškai |
| IT Italiano | LV Latviski |
| RO Română | CS Český |
| NL Nederlands | GA Gaeilge |
| DE Deutsch | HU Magyar |
| SV Svenska | TR Türkçe |
| DA Danish | MT Malti |
| FI Suomi | BG Български |
| HR Hrvatski | RU Русский |
| SL Slovenščina | EL Ελληνικά |

Português

DEFINIÇÃO

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Light-cured masking resin.

INDICATIONS

Masking of dark and stained teeth; metal surfaces during repair of metal-ceramic crowns and other restorations with a metal framework; metal posts and cores in teeth to be restored with metal-free ceramics, indirect composite resin or ceromer; carbon fiber posts; metal surfaces in the dental laboratory.

CONTRAINDICATION

OPAK is contraindicated in patients with history of allergic reaction to any of its components.

COMPOSITION

Bis-GMA, urethane dimethacrylate, catalysts, stabilizers, pigments.

COLORS

B 0.5* - A 3* - Light pink - Dark pink

B and A are basic Vita shades. Vita is a trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

INSTRUCTIONS FOR USE IN THE DENTAL OFFICE

MASKING OF DARK DENTIN

Isolate the operative field (pic. 1) and proceed with: partial removal of dark dentin; application of ANGELUS® 37% PHOSPHORIC ACID for 15 seconds; wash and gently dry; application of bonding agent on dentin and light-curing.

Select an OPAK shade and paint a thin layer on dark dentin and light-cure for 40 seconds (pic. 2/3);

Apply a second layer of OPAK for a more intense masking effect and repeat light-curing. Restore with composite resin (Fig. 4/5).

MASKING OF METAL POSTS AND CORES

Isolate the operative field and increase retentivity of metal (pic. 1/2) or use a primer for metal (MZ PRIMER ANGELUS®);

Apply a bonding agent and light-cure (pic. 3);

Select an OPAK shade and paint a thin layer (pic. 4) and light-cure for 40 seconds;

Apply a second layer of OPAK for a more intense masking effect and repeat light-curing (pic. 5/6);

MASKING OF METAL IN CROWN REPAIRS

Isolate the operative field; bevel the ceramic/resin and roughen the exposed metal with a diamond bur (pic. 1);

Apply 10% ANGELUS® PORCELAIN ETCHANT gel of hydrofluoric acid (HF) on the fracture site of ceramics/composite for 1 minute; wash and dry (pic. 2). **IMPORTANT:** Protect oral tissues against the corrosive and poisonous effects of HF.

Apply a layer of silane (SILANO ANGELUS®) on the ceramic or resin, let it dry for 1 minute and apply the adhesive, then light cure (Fig. 3 a/b).

Select an OPAK shade and paint a thin layer on metal only (pic. 4) and light-cure for 40 seconds;

Apply a second layer of OPAK for a more intense masking effect and repeat light-curing. Restore with composite resin (fig. 5/6).

MASKING OF CARBON FIBER POSTS

After field isolation, dry the cemented post with air (pic. 1);

3. Select the OPAK color, apply a thin layer (Fig. 2), light cure for 40 seconds and perform the filling core (Fig. 3/4).

INSTRUCTIONS FOR USE IN THE DENTAL LAB

Sandblast and apply a metal primer on the metal surface MZ PRIMER ANGELUS® (pics. 1 and 2);

Apply a thin layer of OPAK (light or dark) with a brush (pic. 3) and light-cure for 3 minutes;

Apply extra layers and light-cure each for 3 minutes (pic. 4). Light-cure the final layer for 6 minutes (pic. 5).

PRECAUTIONS

Do not apply thick layers of OPAK;

Isolate area ideally with rubber dam;

Avoid contact of OPAK with eugenol-based substances which inhibit its polymerization;

4. Treat the metal surface prior to application of OPAK (blasting and primer for metal - MZ PRIMER ANGELUS®);

Keep OPAK in refrigerated environment, with temperature lower than 28°C;

Replace syringe cap immediately after use.

Exposure to any source of light may start polymerization;

Contact with non-polymerized OPAK may irritate oral tissues, eyes and skin and cause sensitization in susceptible individuals.

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