Instructions for use

1005102016

Choose your language

- | | |
|-----------------------|-----------------------|
| PT Português | SK Slovensky |
| EN English | PL Polski |
| ES Español | ET Eesti |
| FR Français | LT Lietuviškai |
| IT Italiano | LV Latviski |
| RO Română | CS Český |
| NL Nederlands | GA Gaeilge |
| DE Deutsch | HU Magyar |
| SV Svenska | TR Türkçe |
| DA Danish | MT Malti |
| FI Suomi | BG Български |
| HR Hrvatski | RU Русский |
| SL Slovenščina | EL Ελληνικά |

Português

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English

Adhesion promoter for ceramics and glass fibers

INTRODUCTION

In metal-free dentistry, an adhesive technique is fundamental for the cementation of all-ceramic, indirect composite resin restorations and glass fiber posts.

Adhesion of resin cements to enamel and dentin through acid-etching and application of priming and bonding agents are well known and researched subjects. However, how do resin cements bond to inorganic-based materials such as ceramics and glass fiber posts? Can bonding to cured composite resin be enhanced? The solution is in silanization.

DEFINITION

Silanes are adhesion promoters between organic and inorganic materials. As bi-functional molecules, their silica radicals unite with the silica in the restorative materials and their organic radicals co-polymerize with the organic matrix of the resin cements (methacrylates). Silanes are also called ceramic primers or bonding/coupling agents.

A. CHEMICAL FORMULA

$X-R-Si(OR)_3n$

X: Organofunctional group - provides strong cross-links with resin polymers.

R: Methylene group.

OR: Hydrolizable group - provides strong covalent bonds with hydroxyl groups present in the silica component of ceramics, composites and fiber posts.

Si: Silicon.

n: 0 - 3.

B. PROPERTIES

- Increases wettability of restorative materials in relation to resin cements, enhancing mechanical, physical and chemical adhesion between cement and restoration;
- Increases resistance to dissolution by water at the

cement-restoration interface.

COMPOSITION

- Silane and ethanol.

INDICATIONS

- Internal area of restorations fabricated with ceramic and indirect composite resin for the adhesive cementation technique;
- Glass fiber posts for the adhesive cementation technique;
- Ceramic and indirect composite resin prior to repair with composite resin. If the repair of a restoration is to be made with a fragment of ceramic or indirect composite resin, these should also be silanated.

TECHNIQUE OF USE

A. CEMENTATION OF GLASS FIBER POSTS

1. Clean post with alcohol; dry with air;
2. Apply SILANO ANGELUS® with disposable applicators (APLIK ANGELUS®); wait for 1 (one) minute and gently dry with air (pic. 1);
3. Proceed with cementation following instructions of the adhesive FUSION-DURALINK® ANGELUS® and of the selected resin cement (pic. 2).

B. CEMENTATION OF CERAMIC AND FIBER-REINFORCED INDIRECT COMPOSITE RESIN RESTORATIONS

1. Isolate the margins of the restoration with wax to avoid marginal degradation by the hydrofluoric acid and adhesion of resin cement to the outer surface of the restoration (pic. 3);
2. Apply PORCELAIN CONDITIONING ANGELUS® (10% Hydrofluoric Acid) on the internal area of the restoration; wait for 1 (one) minute* (pic. 4);

* Leucite-based porcelains are etched only for 20 (twenty) seconds.

3. Wash with water and dry with air;
4. Apply SILANO ANGELUS® with disposable applicators (APLIK ANGELUS®); wait for 1 (one) minute and gently dry with air (pic. 5);
5. Proceed with cementation following instructions of the adhesive FUSION-DURALINK® ANGELUS® and of the selected resin cement.

C. REPAIR OF CERAMIC AND INDIRECT COMPOSITE RESIN WITH DIRECT COMPOSITE RESIN

1. Bevel the fractured ceramic surface and roughen the metal (if exposed) or apply a layer of a metal primer (MZ PRIMER ANGELUS®) (pic. 6);
2. Apply PORCELAIN CONDITIONING ANGELUS® (10% Hydrofluoric Acid) on the beveled surface and on the fracture site; wait for 1 (one) minute*, wash abundantly with water and dry (pic. 7);

* Leucite-based porcelains are etched only for 20

(twenty) seconds.

Protect oral tissues against the corrosive and poisonous effects of HF.

3. Apply SILANO ANGELUS® with a brush; wait for 1 (one) minute and gently dry with air (pic. 8);
4. Apply the adhesive FUSION-DURALINK® ANGELUS® bonding agent following its instructions;
5. If there is exposed metal, cover it with a layer of OPAK ANGELUS® masking resin (pic. 9);
6. Restore the fracture site with composite resin (pic. 10).

D. BONDING FRAGMENT OF CERAMIC OR INDIRECT COMPOSITE RESIN

1. Roughen fragment and fracture site;
2. Apply PORCELAIN CONDITIONING ANGELUS® (10% Hydrofluoric Acid) on fracture site and on fragment for 1 (one) minute*; wash and dry (pic. 11);

* Leucite-based porcelains are etched only for 20 (twenty) seconds.

3. Apply SILANO ANGELUS® with disposable applicators (APLIK ANGELUS®); wait for 1 (one) minute and gently dry with air (pic. 12);
4. Cement fragment on restoration following instructions of the adhesive FUSION-DURALINK® ANGELUS® and of the selected resin cement.

WARNINGS | PRECAUTIONS

- Avoid contact with oral tissues, skin and eyes. If accidental contact occurs, wash with water. If necessary, seek medical attention.
- Do not expose near fire or excessive heat. The solution is flammable.
- Keep cap tightly sealed when not in use. The solution is highly volatile.

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