FIBER CAD POST & CORE

ENGLISH

DEFINITION

FIBER CAD Post & Core is a composite of glass fibers and epoxy resin in block formats for 'chairside' use and in disks for laboratory use.

INDICATIONS

Preparation of anatomical and esthetic glass fiber post-and-core by the CAD CAM system.

COMPOSITION

Glass fiber: 75-80%
Epoxy resin: 20-25%

WARNINGS

- Do not use in patients with reported sensitivity to the components of the formula;
- Product of exclusive use by dental technicians and dentists.

ⓘ ***IMPORTANT:** These instructions are based upon clinical and laboratory studies. The success of the procedures depends on the correct diagnosis and careful technique.*

PRECAUTION

- Wear a mask and protective goggles during machining and final adjustments;
- Use rotary instruments (drills); do not use scissors or pliers;
- After the final adjustments and cementation the impermeabilization of the core with a resin flow – by the dentist – is advisable.

ⓘ ***IMPORTANT:** These instructions are based upon clinical and laboratory studies. The success of the procedures depends on the correct diagnosis and careful technique*

STORAGE

Maintain at room temperature (19 to 23 °C) and protected from light.

TECHNICAL INFORMATION

- Density:.....2,15 g/cm³
- Brinell hardness:.....12 to 14 Kg/mm²
- Resistance to flection:.....1100 MPa
- Module of elasticity:.....25 GPa
- Coefficient of thermal expansion (20 - 700 °C):.....10 x 10⁻⁶ °C⁻¹
- Average diameter of fibers:.....33µm

USE TECHNIQUE

LABORATORY

- Scan the preparation of the root in the mouth or of the dental wax/ resin pattern;
- Position FIBER CAD Post & Core in the machining equipment in accordance with the parameters defined by the manufacturer;
- Machine FIBER CAD Post & Core using tungsten carbide drills;
- Remove the post-and-core from the equipment and make the final adjustments with diamond drills.

CLINICAL

A. Adhesive Cementation

Post-and-core preparation

- Place it in the root to check the fit. If adjustments are required, use diamond drills;
- ⓘ ***IMPORTANT:** Do not use scissors or pliers. They mash the fibers and damage the structure;*
- Clean the post-and-core with alcohol and dry it with air;
- Apply a layer of silane and let it to dry for 1 min;
- Apply the chemical adhesive on the radicular part of the structure.

Cementation proper

- Wash the root canal with water; dry with air and absorbent paper points. Rinse irrigating with saline and aspirate, both procedures with Endo Tips Angelus and dry with paper tips;
- Apply the fourth or sixth generation adhesive system (separate flasks) according to the manufacturer's instructions;
- Apply the selected resin cement with a spatula in the root part of the post-and-core and/or in the root canal preparation;
- Place the structure in the root canal, wait for the polymerization and remove the excess of resin cement.

ⓘ ***IMPORTANT:** It is advisable to apply a layer of a flow adhesive resin on the post-and-core for its impermeabilization.*

B. Cementation with zinc phosphate

- Place it in the root to check the fit. If adjustments are required, use diamond points;
- ⓘ ***IMPORTANT:** Do not use scissors or pliers. They mash the fibers and damage the structure.*
- Clean the post-and-core with alcohol and dry it with air;
- Apply the zinc phosphate cement with a spatula in the root part of the post-and-core and/or in the root canal preparation;
- Place the structure in the root canal and keep it under pressure until the final setting time (according to the cement manufacturer);
- Remove the excess.
- ⓘ ***IMPORTANT:** It is advisable to apply a layer of a flow adhesive resin on the post-and-core for its impermeabilization.*

ADVANTAGES OF FIBER CAD POST & CORE

Dentist - Chairside System

- Esthetics: preparation of metal-free prostheses without requiring opacification of the core;
- Strength: withstands mechanical forces near 1000 MPa – 3 times greater than average maximum chewing loads;
- Retention: anatomical fit in root canal promotes high adhesion between post-and-core and root canal;
- Safety: modulus of elasticity similar to dentin´s minimizes risk of

- root fracture;
- Technical Simplicity;
- Same process of root canal preparation, similar to cast metal post technique;
- Easy cementation procedures with adhesive or zinc phosphate cements;
- Easy adjustments if compared to cast metal posts.

Dental Technician - Laboratory System

- Feasibility of CAD-CAM procedure greater than casting process;
- Excellent cost/benefit ratio if compared to cast metal posts;
- Addition of a differentiated and innovative service to the laboratory portfolio.

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