# SILANO

## PORTEGUÊS

### INTRODUÇÃO

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### ENGLISH

#### INTRODUCTION

In metal-free dentistry, an adhesive technique is fundamental for the cementation of all-ceramic, indirect composite resin restorations and glass fiber posts.

Adhesion of resin cements to enamel and dentin through acid-etching and application of priming and bonding agents are well known and researched subjects. However, how do resin cements bond to inorganic-based materials such as ceramics and glass fiber posts? Can bonding to cured composite resin be enhanced? The solution is in silanization.

#### DEFINITION

Silanos are adhesion promoters between organic and inorganic materials. As bi-functional molecules, their silica radicals unite with the silica in the restorative materials and their organic radicals co-polymerize with the organic matrix of the resin cements (methacrylates). Silanes are also called ceramic primers or bonding/coupling agents.

#### CHEMICAL FORMULA

#### X-R-Si(OR)3n

X: Organofunctional group - provides strong cross-links with resin polymers.

R: Methylene group.

OR: Hydrolyzable group - provides strong covalent bonds with hydroxyl groups present in the silica component of ceramics, composites and fiber posts.

Si: Silicon

n: 0 - 3.

#### PROPERTIES

- Increases wettability of restorative materials in relation to resin cements, enhancing mechanical, physical and chemical adhesion between cement and restoration;
- Increases resistance to dissolution by water at the cement-restoration interface.

#### COMPOSITION

- Silane and ethanol.

#### INDICATIONS

- Internal area of restorations fabricated with ceramic and indirect composite resin for the adhesive cementation technique;
- Glass fiber posts for the adhesive cementation technique;
- Ceramic and indirect composite resin prior to repair with composite resin. If the repair of a restoration is to be made with a fragment of ceramic or indirect composite resin, these should also be silanated.

#### TECHNIQUE OF USE

##### A. CEMENTATION OF GLASS FIBER POSTS

- Clean post with alcohol; dry with air;
- Apply ANGELUS® SILANO with disposable applicators (ANGELUS® APLIK); wait for 1 (one) minute and gently dry with air (pic. 1);
- Proceed with cementation following instructions of the adhesive ANGELUS® FUSION-DURALINK® and of the selected resin cement (pic. 2).

##### B. CEMENTATION OF CERAMIC AND FIBER-REINFORCED INDIRECT COMPOSITE RESIN RESTORATIONS

- Isolate the margins of the restoration with wax to avoid marginal degradation by the hydrofluoric acid and adhesion of resin cement to the outer surface of the restoration (pic. 3);
- Apply ANGELUS® PORCELAIN ETCHANT (10% Hydrofluoric Acid) on the internal area of the restoration; wait for 1 (one) minute\* (pic. 4);
- \* Leucite-based porcelains are etched only for 20 (twenty) seconds.

3. Wash with water and dry with air;

4. Apply ANGELUS® SILANO with disposable applicators (ANGELUS® APLIK); wait for 1 (one) minute and gently dry with air (pic. 5);

5. Proceed with cementation following instructions of the adhesive ANGELUS® FUSION-DURALINK® and of the selected resin cement.

##### C. REPAIR OF CERAMIC AND INDIRECT COMPOSITE RESIN WITH DIRECT COMPOSITE RESIN

1. Bevel the fractured ceramic surface and roughen the metal (if exposed) or apply a layer of a metal primer (ANGELUS® MZ PRIMER) (pic. 6);
2. Apply ANGELUS® PORCELAIN ETCHANT (10% Hydrofluoric Acid) on the beveled surface and on the fracture site; wait for 1 (one) minute\*, wash

abundantly with water and dry (pic. 7);

\* Leucite-based porcelains are etched only for 20 (twenty) seconds.

④ Protect oral tissues against the corrosive and poisonous effects of HF.

3. Apply ANGELUS® SILANO with a brush; wait for 1 (one) minute and gently dry with air (pic. 8);

4. Apply the adhesive ANGELUS® FUSION-DURALINK® bonding agent following its instructions;

5. If there is exposed metal, cover it with a layer of ANGELUS® OPAK masking resin (pic. 9);

6. Restore the fracture site with composite resin (pic. 10).

##### D. BONDING FRAGMENT OF CERAMIC OR INDIRECT COMPOSITE RESIN

1. Roughen fragment and fracture site;
2. Apply ANGELUS® PORCELAIN ETCHANT (10% Hydrofluoric Acid) on fracture site and on fragment for 1 (one) minute\*, wash and dry (pic. 11);

\* Leucite-based porcelains are etched only for 20 (twenty) seconds.

3. Apply ANGELUS® SILANO with disposable applicators (ANGELUS® APLIK); wait for 1 (one) minute and gently dry with air (pic. 12);

4. Cement fragment on restoration following instructions of the adhesive ANGELUS® FUSION-DURALINK® and of the selected resin cement.

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- Avoid contact with oral tissues, skin and eyes. If accidental contact occurs, wash with water. If necessary, seek medical attention.
- Do not expose near fire or excessive heat. The solution is flammable.
- Keep cap tightly sealed when not in use. The solution is highly volatile.

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