PORTUGUÊS

DEFINIÇÃO

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ENGLISH

DEFINITION

System of fiberglass intraradicular pin.

INDICATIONS

Support of prosthetic crowns and restorations.

COMPOSITION

Pins and sleeve (% in weight): Fiberglass.....80%
Epoxy resin.....20%

Burs.....Stainless steel

TECHNIQUE OF USE

A. PREPARATION OF ROOT CANAL

1. Remove the filling of the canal partially with heated instruments to the pre-determined depth. Leave at least 3 to 5 mm of filling in the root apex;
2. Prepare the duct with the Splendor bur;
3. X-ray to check the preparation of the canal and condition of the filling.

B. PREPARE THE SPLENDOR SYSTEM FOR CEMENTING

1. Position the pin in the duct until the end of the preparation;
2. Insert the sleeve on the pin and position it as near the apex as possible with slight pressure;
3. Check the occlusion and cut the pin and sleeve unit at the height defined;

4. Remove the pin and sleeve unit from the duct;
5. Clean the sleeve and pin with alcohol and dry them with air jets;
6. Apply a layer of silane (SILANO ANGELUS®) and leave to dry for 1 minute;
7. Apply the adhesive system on the pin and sleeve in accordance with the instructions for use of the adhesive used.

C. CEMENTING

1. Wash the conduct with water jets and dry it with absorbing paper tips;
 2. Prepare the duct in accordance with the instructions for use of the cement used;
 3. Spatulate the cement and fill the duct;
 4. Position initially the pin in the duct and then the sleeve on it;
 5. Remove the excesses of cement and wait for its polymerization.
- IMPORTANT:** Cements of zinc phosphate and glass ionomer (conventional and modified by resin) can also be used, but their mechanical properties are inferior to those of resinous cements for cementing. The self-adhesive cements dispense with the steps of acid conditioning and application of adhesive.

D. PREPARATION OF THE CORONARY PART OF THE CORE

1. Construct the coronary part of the core with composite resin.

WARNINGS | PRECAUTIONS

- Use diamond tips at high rotation and under cooling to cut the fiberglass pins;
- Do not use scissors or pliers to cut the fiberglass pins;
- Wear safety googles and a mask to avoid inhaling particles generated during cutting;
- Ensure that the fiberglass pins were totally covered by the cementing agent or by the composite resin. Exposed fibers absorb water, causing alterations to the mechanical properties of the pin.

STORAGE

- Keep in a dry, ventilated place out of the sunlight.

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10506250 - 0721082023

For the digital instructions,
in other languages please access this QR code.

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ATTENTION: This product must be used according to the instructions in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use.

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Polypropylene.

INSTRUCTIONS FOR USE

1. Disinfect the tips by immersion in 1% sodium hypochlorite or 70% ethanol for 30 minutes;
2. Wash with water to remove debris and dry with air;
3. Attach tips to cannulas for suction.

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