# OPAK

## PORTUGUÊS

### INTRODUÇÃO

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## ENGLISH

### INTRODUCTION

With the aesthetic needs improvement on Dentistry, opakers composites were developed to cover metal and dark surfaces.

### DEFINITION

OPAK is a light-cured masking resin for the covering of dark dentin, metal and carbon fiber posts. And for the covering of metal surfaces in the dental laboratory.

### INDICATIONS

- a. Dark and stained teeth.
- b. Metal surfaces during repair of metal-ceramic crowns and other restorations with a metal framework.
- c. Metal posts and cores in teeth to be restored with metal-free ceramics, indirect composite resin or ceromer.
- d. Carbon fiber posts
- e. Metal surfaces in the dental laboratory.

### CHARACTERISTICS

- a. High masking capacity.
- b. Low viscosity: can be applied with brushes.
- c. Excellent adhesion.
- d. Compatible with all bonding agents and composite resins.

### COMPOSITION

Bis-GMA, urethane dimethacrylate, catalysts, stabilizers, pigments.

### COLORS

B 0.5*	A 3*	Light	Dark
--------	------	-------	------

\*B and A are basic Vita shades. Vita is a trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

### PRECAUTIONS

- a. Do not apply thick layers of OPAK.
- b. Isolate area ideally with rubber dam.
- c. Avoid contact of OPAK with eugenol-based substances which inhibit its polymerization.
- d. Prepare the metal surface before the application of OPAK (sandblasting and metal primer - MZ PRIMER ANGELUS®).
- e. Keep OPAK in refrigerated environment, with temperature lower than 28°C.
- f. Replace syringe cap immediately after use. Prolonged exposure to any source of light may start polymerization.
- g. Contact with non-polymerized OPAK may irritate oral tissues, eyes and skin and cause sensitization in susceptible individuals.

### WARNINGS

- a. Replace syringe cap immediately after use. Prolonged exposure to any source of light may start polymerization.

b.Contact with non-polymerized OPAK may irritate oral tissues, eyes and skin and cause sensitization in susceptible individuals.

c.Keep OPAK in temperatures between 2°C and 28°C.

## **CONTRAINDICATION**

OPAK is contraindicated in patients with history of allergic reaction to any of its components.

## **INSTRUCTIONS FOR USE IN THE DENTAL OFFICE**

### **A.MASKING OF DARK DENTIN**

1.Isolate the operative field (pic. 1).

2.Proceed with:

- Partial removal of dark dentin.
- Application of ANGELUS® 37% PHOSPHORIC ACID for 15 seconds; wash and gently dry.
- Application of bonding agent on dentin and light-curing.

3.Select an OPAK shade (B 0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on dark dentin (pic. 2).

4.Light-cure for 40 seconds (pic. 3).

5.If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B 0.5). Repeat light-curing (pic. 4).

6.Apply a composite resin and finish the restoration (pic. 5).

### **B.MASKING OF METAL POSTS AND CORES**

1.Isolate the operative field (pic. 1).

2.Increase retentivity of metal (roughening with burs or sandblasting) or use a primer for metal (MZ PRIMER ANGELUS®) (pic. 2).

3.Apply a bonding agent and light-cure (pic. 3).

4.Select an OPAK shade (B 0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on metal (pic. 4).

5.Light-cure for 40 seconds.

6.If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B 0.5). Repeat light-curing (pic. 5).

7.Final result (pic. 6).

### **C. MASKING OF METAL IN CROWN REPAIRS**

1.Isolate the operative field; bevel the fractured ceramic or indirect resin surface and roughen the exposed metal with a diamond bur (pic. 1).

2.Apply 10% ANGELUS® PORCELAIN ETCHANT gel of hydrofluoric acid (HF) on the beveled surface and on the fracture site of ceramics or composite for 1 minute; wash abundantly with water and dry (pic. 2).

*IMPORTANT: Protect oral tissues against the corrosive and poisonous effects of HF.*

3.Apply SILANE ANGELUS® with brush on ceramics or composite; wait for 1 minute and gently dry with air (pic. 3a).

4.Apply a bonding agent and light-cure (pic. 3b).

5.Select an OPAK shade (B0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on metal only. Light-cure for 40 seconds (pic. 4).

6.If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B0.5). Repeat light-curing (pic. 5).

7.Apply a composite resin and finish the repair (pic. 6).

#### **D.MASKING OF CARBON FIBER POSTS**

- 1.After field isolation, dry the cemented post with air (pic. 1).
- 2.Select an OPAK shade and paint the post with a thin layer (0.3 mm) (pic. 2).
- 3.Light-cure all surfaces for 40 seconds each. If a more intense masking effect is needed, apply a second layer of OPAK. Repeat light-curing (pic. 3).
- 4.Composite resin core build-up ready (pic. 4).

#### **INSTRUCTIONS FOR USE IN THE DENTAL LAB**

- 1.Sandblast and apply a metal primer on the metal surface MZ PRIMER ANGELUS® (pics. 1 and 2).
- 2.Apply a thin layer of the selected OPAK (light or dark) with a brush (pic. 3).
- 3.Light-cure for 3 minutes.
- 4.Apply two extra layers; light-cure each for 3 minutes (pic. 4).
- 5.Light-cure the final layer for 6 minutes (pic. 5).

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