
MTA-Fillapex

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35. GLASSMANN, G. Bioactive endodontic obturation: Combining the new with the tried and true MTA Fillapex and Continuous Wave of Condensation; ROOTS, 3; 2013
36. GOMES FILHO, J E; Rat tissue reaction to MTA-Fillapex Dental Traumatology 2011; doi: 10.1111/j.1600-9657.2011.
37. GOMES FILHO, J. E.; et al. Effect of MTA-based sealer on the healing of periapical lesions. J Appl Oral Sci.2013;21(3):235-42
38. GOMES-FILHO J.E, MOREIRA J.V, WATANABE S., LODI C.S., CINTRA L.T.S, DEZAN JUNIOR, E., BERNABÉ P.F.E, NERY M.J, OTOBONI FILHO, J.A Sealability of MTA and calcium hydroxide containing Sealers. J Appl Oral Sci. 2012;20(3):347-51
39. GOMES-FILHO JE, WATANABE S, BERNABÉ PFE, COSTA MTM. A mineral trioxide aggregate sealer stimulated mineralization. Journal of Endodontics 35:256- 60 (2009)
40. GÜMRÜ, B.; ET AL. Evaluation of The Radiopacity of a MTA-Based Root-Canal Filling Material using Digital Radiography. Journal of Marmara University Institute of Health Sciences Volume: 3, Number: 1, 2013
41. IOANNIDIS, K.; ET AL. Spectrophotometric analysis of crown discoloration induced by MTA- and ZnOE-based sealers. J. Appl. Oral Sci. vol.21 no.2 Bauru Mar./Apr. 2013
42. KESKIN, S.; ET AL. Efficacy of Different Activation Regimes of Chitosan. International Dental Journal © 2013 FDI World Dental Federation
43. KUGA , M .C.; ET AL. Effects of calcium hydroxide addition on the physical and chemical properties of a calcium silicate-based sealer. J Appl Oral Sci. 2014 May-Jun; 22(3): 180–184
44. KUGA, M. C.; ET AL. Evaluation of the pH, calcium release and antibacterial activity of MTA Fillapex. Rev Odontol UNESP. 2013 Sept-Oct; 42(5): 330-335
45. KUGA,M.C. et al. Hydrogen ion and calcium releasing of MTA-Fillapex® and MTA-based formulations. RSBO. 2011 Jul-Sep;8(3):271-6
46. MACHADO, G.L.N. Estudo Comparativo Da Infiltração Apical Entre Um Cimento A Base De MTA E Um Cimento Resinoso. Monografia apresentada ao Programa de Especialização em Endodontia do ICS – FUN Ilhéus-BA, como parte dos requisitos para obtenção do título de Especialista. Ilhéus – BA, 2013
47. MADANI, Z .S.; ET AL. Comparative evaluation of antimicrobial activity of two root canal sealers: MTA Fillapex and AH26. Minerva Stomatol, 2014. 63; 1-2
48. MALKA, V.B. Radiopacidade de cimentos endodonticos: comparação entre dois métodos in vitro.Trabalho de conclusão de curso (graduação). Universidade Federal do Rio Grande do Sul. Porto Alegre 2012
49. MARCIANO, M.A. Physical And Chemical Properties Of Experimental Portland-based Sealers Saturday, March 23, 2013: 2 p.m. - 3:15 p.m. Location: Hall 4 (Washington State Convention Center)IADR, MARCH, 2013

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Sealability of MTA and calcium hydroxide-containing sealers

João Eduardo GOMES-FILHO², Jaqueline Viana MOREIRA⁵, Simone WATANABE⁴, Carolina Simonetti LODI⁴, Luciano Tavares Angelo CINTRA³, Eloi DEZAN JUNIOR², Pedro Felício Estrada BERNABÉ¹, Mauro Juvenal NERY², José Arlindo OTOBONI FILHO²

1- DDS, MSc, PhD, Full Professor, Department of Endodontics, Araçatuba School of Dentistry, UNESP - Univ. Estadual Paulista, Araçatuba, SP, Brazil.

2- DDS, MSc, PhD, Associate Professor, Department of Endodontics, Araçatuba School of Dentistry, UNESP - Univ. Estadual Paulista, Araçatuba, SP, Brazil.

3- DDS, MSc, PhD, Assistant Professor, Department of Endodontics, Araçatuba School of Dentistry, UNESP - Univ. Estadual Paulista, Araçatuba, SP, Brazil.

4- DDS, MSc, Graduate Program in Pediatrics, Araçatuba School of Dentistry, UNESP - Univ. Estadual Paulista, Araçatuba, SP, Brazil.

5- Undergraduate student, Araçatuba School of Dentistry, UNESP - Univ. Estadual Paulista, Araçatuba, SP, Brazil.

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Received: October 2, 2010 - Modification: August 16, 2011 - Accepted: September 1, 2011

ABSTRACT

Objectives: The aim of this study was to evaluate the apical sealability of Fillapex[®], Endo-CPM-Sealer[®] and Sealapex[®]. **Material and Methods:** Ninety-four freshly extracted single-rooted teeth were selected and decoronated. All teeth were radiographed to confirm the existence of a single and straight root canal, which was prepared using Protaper Universal and 2.5% sodium hypochlorite. The teeth were randomly divided in groups of 10 specimens each according to the sealer, and the canals were filled using the single cone technique and one of the sealers. Four additional teeth were used as controls. The teeth were submitted to dye leakage with Rhodamine B for 24 h but using vacuum on the initial 15 min. Thereafter, they were cut longitudinally and the leakage was measured in a linear fashion from apex to crown. Data were analyzed by ANOVA and Tukey's tests at 5% significance level. **Results:** Fillapex[®] and Sealapex[®] showed significantly less dye leakage than Endo-CPM-Sealer[®] ($p < 0.05$). **Conclusions:** It was concluded that Fillapex[®] and Sealapex[®] were able to prevent apical dye leakage differently from Endo-CPM-Sealer[®].

Key words: Root canal filling material. Leakage. Fillapex.

Original Research Article

Hydrogen ion and calcium releasing of MTA Fillapex[®] and MTA-based formulations

Milton Carlos Kuga¹
Edson Alves de Campos¹
Paloma Hernandez Viscardi²
Paula Zapparoli Carrilho²
Fernanda Castilho Xavier²
Nayara Pereira Silvestre²

Corresponding author:

Milton Carlos Kuga
Avenida Saul Silveira, n. 5-01.
CEP 17018-260 – Bauru – SP – Brasil
E-mail: kuga@foar.unesp.br

¹ School of Dentistry of Araraquara, Sao Paulo State University – Araraquara – SP – Brazil.

² School of Dentistry, University Center of Rio Preto – Sao Jose do Rio Preto – SP – Brazil.

Received for publication: January 12, 2011. Accepted for publication: February 23, 2011.

Keywords:

MTA; pH; calcium; sealer.

Abstract

Introduction: MTA is composed of various metal oxides, calcium oxide and bismuth. It has good biological properties and is indicated in cases of endodontic complications. Several commercial formulations are available and further studies are necessary to evaluate these materials. **Objective:** To evaluate pH and calcium releasing of MTA Fillapex[®] compared with gray and white MTA. **Material and methods:** Gray and white MTA (Angelus) and MTA Fillapex[®] (Angelus) were manipulated and placed into polyethylene tubes and immersed in distilled water. The pH of these solutions was measured at 24 hours, 7 days and 14 days. Simultaneously, at these same aforementioned periods, these materials' calcium releasing was quantified, through atomic absorption spectrophotometry. The results were submitted to ANOVA, with level of significance at 5%. **Results:** Concerning to pH, the materials present similar behaviors among each other at 24 hours ($p > 0.05$). At 7 and 14 days, MTA Fillapex[®] provided significantly lower pH values than the other materials ($p < 0.05$). Regarding to calcium releasing, at 24 hours and 7 days, MTA Fillapex[®] provided lower releasing than the other materials ($p < 0.05$). After 14 days, differences were found between MTA Fillapex[®] and gray MTA ($p < 0.05$). **Conclusion:** All materials showed alkaline pH and calcium releasing, with significantly lower values for MTA Fillapex[®] sealer.

ISSN:
Printed version: 1806-7727
Electronic version: 1984-5685
RSBO. 2011 Oct-Dec;8(4):404-9

Original Research Article

Comparison of MTA Fillapex radiopacity with five root canal sealers

Ana Paula Meirelles Vidotto¹
Rodrigo Sanches Cunha²
Eduardo Gregatto Zeferino¹
Daniel Guimarães Pedro Rocha¹
Alexandre Sigrist de Martin¹
Carlos Eduardo da Silveira Bueno¹

Corresponding author:

Ana Paula Meirelles Vidotto
Avenida Santa Isabel, n.º 241
CEP 13084-755 – Campinas – SP – Brasil
E-mail: ana_vidotto@hotmail.com

¹ São Leopoldo Mandic Research Center – Campinas – SP – Brazil.

² University of Manitoba – Winnipeg – MB – Canada.

Received for publication: February 4, 2011. Accepted for publication: March 30, 2011.

Keywords: obturation;
root canal sealer;
radiopacity.

Abstract

Introduction: The endodontic sealer is a filling material whose physicochemical properties are mandatory for the achievement of endodontic therapy final goal. An ideal endodontic sealer should have some properties, including radiopacity. **Objective:** This study compared MTA Fillapex™ radiopacity with the radiopacity of five others endodontic sealers: Endométhasone-N™, AH Plus™, Acroseal™, Epiphany SE™ and RoekoSeal™. **Material and methods:** Five cylindrical samples of each sealer were used, constructed with the aid of a matrix. On an occlusal film, a sample of each sealer was placed along with an aluminum stepwedge and five radiographic shots were taken. The radiographic images were digitized and each sample's gray scales were compared with each shade of the aluminum stepwedge, by using software. **Results:** The results, in decreasing order of radiopacity, were: AH Plus™ was statistically the most radiopaque sealer (9.4 mm Al), followed by Epiphany SE™ (7.8 mm Al), MTA Fillapex™ (6.5 mm Al), RoekoSeal™ (5.8 mm Al), Endométhasone-N™ (4.5 mm Al), and Acroseal™, the least statistically radiopaque (3.5 mm Al). **Conclusion:** It can be concluded that MTA Fillapex™ was the third most radiopaque sealer among all tested sealers. Also, MTA Fillapex™ has the radiopacity degree in agreement with ADA specification No. 57 (1983).

Dental Traumatology

Dental Traumatology 2011; doi: 10.1111/j.1600-9657.2011.01096.x

Rat tissue reaction to MTA FILLAPEX[®]

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Department of Endodontics, Araçatuba School of Dentistry, University of Estadual Paulista, São Paulo, Brazil

Correspondence to: Dr João Eduardo Gomes-Filho, Department of Endodontics, Araçatuba School of Dentistry, University of Estadual Paulista, R. José Bonifácio, 1193, Araçatuba, São Paulo, Brazil
Tel.: +0055 18 36363252
Fax: +0055 18 36363279
e-mail: joao@foa.unesp.br

Accepted 6 November, 2011

Abstract – The aim of this study was to evaluate the rat subcutaneous tissue reaction to implanted polyethylene tubes filled with mineral trioxide aggregate (MTA) FILLAPEX[®] compared to the reaction to tubes filled with Sealapex[®] or Angelus MTA[®]. These materials were placed in polyethylene tubes and implanted into the dorsal connective tissue of Wistar rats for 7, 15, 30, 60, and 90 days. The specimens were stained with hematoxylin and eosin or Von Kossa or left unstained for examination under polarized light. Qualitative and quantitative evaluations of the reaction were performed. All materials caused moderate reactions after 7 days, which decreased with time. The reactions were moderate and similar to that evoked by the control and Sealapex[®] on the 15th day. MTA FILLAPEX[®] and Angelus MTA caused mild reactions beginning after 15 days. Mineralization and granulation birefringent to polarized light were observed with all materials. It was concluded that MTA FILLAPEX[®] was biocompatible and stimulated mineralization.

Mineral Trioxide Aggregate–based Endodontic Sealer Stimulates Hydroxyapatite Nucleation in Human Osteoblast-like Cell Culture

Loise Pedrosa Salles, MSc,^{*†} Ana Livia Gomes-Cornélio, MSc,^{*} Felipe Coutinho Guimarães, MSc,[†] Bruno Schneider Herrera, PhD,[‡] Sonia Nair Bao, PhD,[†] Carlos Rossa-Junior, PhD,[‡] Juliane Maria Guerreiro-Tanomaru, PhD,^{*} and Mario Tanomaru-Filho, PhD^{*}

Abstract

Introduction: The main purpose of this study was to evaluate the biocompatibility and bioactivity of a new mineral trioxide aggregate (MTA)-based endodontic sealer, MTA Fillapex (MTA-F; Angelus, Londrina, Brazil), in human cell culture. **Methods:** Human osteoblast-like cells (Saos-2) were exposed for 1, 2, 3, and 7 days to MTA-F, Epiphany SE (EP-SE; SybronEndo, Orange, CA), and zinc oxide–eugenol sealer (ZOE). Unexposed cultures were the control group (CT). The viability of the cells was assessed by MTT assay and the morphology by scanning electron microscopy (SEM). The bioactivity of MTA-F was evaluated by alkaline phosphatase activity (ALP) and the detection of calcium deposits in the culture with alizarin red stain (ARS). Energy-dispersive X-ray spectroscopy (EDS) was used to chemically characterize the hydroxyapatite crystallites (HAP). Saos-2 cells were cultured for 21 days for ARS and SEM/EDS. ARS results were expressed as the number of stained nodules per area. Statistical analysis was performed with analysis of variance and Bonferroni tests ($P < .01$). **Results:** MTA-F exposure for 1, 2, and 3 days resulted in increased cytotoxicity. In contrast, viability increased after 7 days of exposure to MTA-F. Exposure to EP-SE and ZOE was cytotoxic at all time points. At day 7, ALP activity increase was significant in the MTA-F group. MTA-F presented the highest percentage of ARS-stained nodules (MTA-F > CT > EP-SE > ZOE). SEM/EDS analysis showed hydroxyapatite crystals only in the MTA-F and CT groups. In the MTA-F group, crystallite morphology and chemical composition were different from CT. **Conclusions:** After setting, the cytotoxicity of MTA-F decreases and the sealer presents suitable bioactivity to stimulate HAP crystal nucleation. (*J Endod* 2012; ■:1–6)

Key Words

Bioactivity, biocompatibility, hydroxyapatite, mineral trioxide aggregate sealer

Mineral trioxide aggregate (MTA) emerged as the material of choice for root perforation repairs and root-end fillings in the 90s, a revolutionary period marked by many advances in endodontics (1). MTA was developed at Loma Linda University and received approval from the Food and Drug Administration for human use in 1998 (2, 3). Since then, MTA has shown excellent biological properties in several *in vivo* and *in vitro* studies (4–9). In cell culture systems, for example, MTA has been shown to enhance proliferation of periodontal ligament fibroblasts (6), to induce differentiation of osteoblasts (7, 8), and to stimulate mineralization of dental pulp cells (9). This biocompatibility and bioactive potential raised the interest of scientists worldwide to improve the handling characteristics and some physicochemical properties of MTA with the intention of expanding its applicability in endodontics. Consequently, new MTA-based root-end filling cements and root canal sealers have been proposed (10–12), such as MTA Fillapex (MTA-F; Angelus, Londrina, Brazil).

The new MTA-based sealers reflect a current requirement to have materials for endodontic therapy that are able to stimulate the healing process of periapical tissues, instead of merely biocompatible or inert materials. As a result, MTA-F represents the effort in combining a material of excellent biological properties as MTA with resins and other components to improve diverse required properties of an endodontic sealer including adhesiveness, dimensional stability, working time, radiopacity, flow, and antibacterial effects. According to the manufacturer's information, MTA-F is composed of salicylate resin, resin diluent, natural resin, bismuth oxide as radiopacifying agent, silica nanoparticles, MTA, and pigments. The MTA itself consists of fine hydrophilic particles of tricalcium silicate, tricalcium aluminum oxide, tricalcium oxide, gypsum (calcium sulfate dihydrate), and other mineral oxides (3). Gypsum is an important determinant of setting time. MTA cements generally contain less gypsum to allow more handling time. Unfortunately, MTA-F data sheet lacks details about the natural resin, pigments, and diluents composition.

It is important to investigate if the combination of these resins and other constituents influence the bioactive potential of MTA in the new endodontic sealer. Therefore, the main purpose of this study was to evaluate the biocompatibility and the bioactivity of MTA-F in stimulating mineralization in Saos-2 cell culture compared with Epiphany SE

[Start](#) | [Browse by Day](#) | [Author Index](#) | [Keyword Index](#)

721 Evaluation of a MTA-based sealer histocompatibility in subcutaneous tissue

Thursday, June 21, 2012: 11:45 a.m. - 1 p.m.

Location: Poster Hall (Convention Center)

Presentation Type: Poster Session

R.V. FERINO, M. TANOMARU-FILHO, G.F. SILVA, E. SASSO-CERRI, J.M. GUERREIRO-TANOMARU, and P.S. CERRI, *Dentística Restauradora, FACULDADE DE ODONTOLOGIA DE ARARAQUARA, UNESP Univ Estadual Paulista, Câmpus Araraquara, Araraquara / SP, Brazil*

Objective: Mineral Trioxide Aggregate (MTA) has been widely used in endodontic treatment due to its good marginal adaptation and biocompatibility. In attempt to associate physicochemical properties of a root canal sealer with MTA biocompatibility, some modifications have been made with the aim to improve its use as a root canal sealer. Thus, the tissue reaction promoted by a MTA-based sealer (MTA-Fillapex) in rat subcutaneous was investigated by morphological and morphometric analyses

Method: Eighty rats were distributed into 4 groups (n=20); in each animal, a polyethylene tube filled with MTA-Fillapex, MTA, AH-Plus or Fill Canal was implanted in the dorsal subcutaneous. After 7, 15, 30 and 60 days, the tubes surrounded by connective tissue were removed, fixed and embedded in paraffin. In the HE-stained sections, the numerical density of inflammatory cells (IC) in the capsule was evaluated and statistical analyses were performed using ANOVA and Tukey's test ($p \leq 0.05$). Sections were also submitted to von Kossa method for detection of calcified structures

Result: At 7 and 15 days, the number of IC was significantly higher in the capsule of the MTA-Fillapex in comparison to other materials. However, significant reduction in the number of IC was verified in the capsule of the MTA-Fillapex in the 30 and 60 days when compared to the initial periods. At 60 days, the inflammatory reaction promoted by MTA-Fillapex was similar to MTA and significantly lower than Fill Canal. Otherwise, a gradual and significant increase was observed in the number of IC in the capsule adjacent to the Fill Canal. von Kossa-positive structures were observed in the capsule adjacent to the MTA-Fillapex, MTA and AH-Plus; positive structures were not seen in the capsule surrounding the Fill Canal

Conclusion: The results indicate that the biocompatibility of MTA-Fillapex is similar to MTA

Keywords: Biocompatibility, Dental materials, Endodontics and Root canal fillings

See more of: [Biocompatibility and Biologic Effects II](#)
See more of: [Dental Materials 5: Biocompatibility and Biologic Effects](#)

[<< Previous Abstract](#) | [Next Abstract](#)

[Abstracts](#) [Author Index](#) [Keyword Index](#) [Search](#) [Help](#) [Main Menu](#)

89TH GENERAL SESSION & EXHIBITION OF THE IADR • SAN DIEGO CONVENTION CENTER • SAN DIEGO, CALIFORNIA, USA

1164 MTA FILLAPEX[®] STIMULATED MINERALIZATION

Thursday, March 17, 2011: 3:30 p.m. - 4:45 p.m.

Location: Hall C (San Diego Convention Center)

Presentation Type: Poster Session

J.E. GOMES-FILHO¹, S. WATANABE¹, L.T.A. CINTRA¹, C.S. LODI², E. DEZAN, Jr.³, M.J. NERY⁴, and P.F.E. BERNABE^{1, 2}, *¹Endodontics, Universidade Est. Paulista Julio Mesquita, Araçatuba, Brazil, ²Pediatric Dentistry, Universidade Est. Paulista Julio Mesquita, Araçatuba, Brazil, ³Endodontics, Universidade Est. Paulista Julio Mesquita, Araçatuba, Brazil, ⁴Endodontics, Universidade Est. Paulista Julio Mesquita, Araçatuba, Brazil*

Objectives: The aim of this study was to evaluate the rat subcutaneous tissue response to implanted polyethylene tubes filled with MTA Fillapex[®] (Angelus, Londrina, Brazil) compared Endo-CPM-Sealer[®] (Portland Cement Modified Sealer) (EGEO S.R.L., Buenos Aires, Argentina) and Sealapex[®] (SybronEndo, Glendora, CA). **Methods:** These materials were placed in polyethylene tubes and implanted into dorsal connective tissue of Wistar rats for 7, 15, 30, 60, and 90 days. The specimens were prepared to be stained with hematoxylin and eosin or Von Kossa or not stained for polarized light. Qualitative and quantitative evaluations of the reaction were performed. **Results:** All materials caused moderate reactions at 7 days that decreased with time. The response was similar to the control on the 15th day with Endo-CPM-Sealer[®] and Sealapex[®]. MTA Fillapex[®] caused mild reactions from 15 days on. Mineralization and granulations birefringent to the polarized light were observed with all materials. **Conclusions:** It was possible to conclude that MTA Fillapex[®] was biocompatible and stimulated mineralization. Supported by FAPESP #2009/11194-5 and #2009/00380-2

Keywords: Animal, Endodontics, MTA, Root canal fillings and Wound healing

See more of: [Dental Materials: Biocompatibility and Biologic Effects II](#)
See more of: [Dental Materials 5: Biocompatibility and Biologic Effects](#)

[<< Previous Abstract](#) | [Next Abstract >>](#)

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☎ 43 9 9919 0244

0800 727 3201

+55 43 2101 3200

 /angelusodonto
/angelusodontoglobal

 /AngelusOdonto
/AngelusOdontoGlobal

www.angelus.ind.br

 angelus[®]

 angelus[®]
Japan

 angelus[®]
USA

PRIMA DENTAL
by  angelus[®]