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PORTUGUÊS

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ENGLISH

THUMB SUCKING

CAUSES: Hunger, insecurity; desire to attract attention motivated by the birth of a sibling, or jealousy arising by any change in the home, such as parental separation, family conflict; shorter duration of breastfeeding.

IMPACTS: Compromising of occlusion (anterior open bite), space between teeth, tightening of the jaw, cross bite, projection of upper incisors, changes in psychological development, speech changes; breathing through mouth, bruxism. Can also compromise some factors such as aesthetics, self-esteem, proper nutrition and even systemic balance.

TREATMENT: All approaches must be made with understanding and patience. Stances that lead patients to feel guilty for having such a habit should not be taken as this can cause trauma to the child. Monitoring must be carried out by the dentist in conjunction with an orthodontist, speech therapist, and sometimes with the help of an ENT specialist and a psychologist. The act of thumb sucking has a multifactorial etiology, and the success of its treatment depends on understanding and cooperation between professionals, family members and the patient. The approach should also include the treatment of associated disorders.

HABIT-BREAKER is an innovative product, easy to use, that acts as a tool in psychological treatment of the child both professionally and by parents. The understanding and cooperation of the child are key to changing habits, so attractive, playful designs will serve as an incentive and reinforcement for successful treatment. HABIT-BREAKER should only be applied to the fingers used for sucking to increase the child's focus of attention. It may be changed periodically with every step forward in breaking the habit. HABIT-BREAKER works in cases of unconscious thumb sucking which tends to become an increasingly frequent and intensive habit. In these cases, making the child conscious of the act of putting fingers in the mouth can be one of the first and most important steps for the psychological treatment of the child. HABIT-BREAKER is an aid in this process, but should not be considered as a means of definitive treatment. The therapeutic approach should always be based on etiology, and on the involvement of the patient and family.

APPLICATION: 1. The application must be performed by a Dental Surgeon or a parent or guardian, so it will not be seen as a play, but as a treatment; 2. The figure should be selected, cut and the film covering the tattoo must be removed; 3. The finger where the HABIT-BREAKER sticker is to be applied must be clean and the figure should be positioned, pressing hard on the print against the skin; 4. Wipe the paper lightly with a damp or cotton wool, still pressing down. Stop when the paper has turned clear; 5. Remove the damp paper carefully. It should come

off without any effort. If it is difficult to remove, it is because it wasn't sufficiently damp. If this happens, start wiping again with the cotton wool; 6. Leave the HABIT-BREAKER to dry for a few minutes. It can last for up to five days; 7. To remove it, wipe it with cotton wool dampened with a gentle make-up remover.

WARNING - IMPORTANT: The HABIT-BREAKERS may cause allergic reactions. If this occurs, discontinue use immediately. For sensitive skin, test on the back of your hand first. In the event of irritation, discontinue use.

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