

# CONTACTO PAPEL

## Instructions for use

### Choose your language

- ENGLISH - CONTACTO PAPEL
- ESPAÑOL - CONTACTO PAPEL
- PORTUGUÊS - CONTACTO PAPEL
- TÜRKÇE - CONTACTO PAPEL



## ENGLISH - CONTACTO PAPEL

### COMPOSITION:

Special paper, vegetal wax and oils, pigments, paraffin and petroleum jelly.

### TECHNIQUE OF USE:

Place articulating paper on the dry occlusal surface and articulate to obtain the occlusal marks.

### PRECAUTION:

Store in a cool dry.

## ESPAÑOL - CONTACTO PAPEL

### COMPOSICIÓN:

Papeles especiales, ceras y aceites vegetales, pigmentos, parafina, vaselina.

### TÉCNICA DE USO:

Posicione el papel articular sobre la superficie oclusal seca y articule para marcar los puntos de contacto.

### PRECAUCIÓN:

Almacenar en un local seco y aireado.

## PORTUGUÊS - CONTACTO PAPEL

### COMPOSIÇÃO:

Papéis especiais, ceras e óleos vegetais, pigmentos, parafina e vaselina.

### TÉCNICA DE USO:

Posicione o carbono na superfície oclusal seca e articule para obter o registro dos pontos de contato.

### PRECAUÇÃO:

Armazenar em local seco e arejado.



# TÜRKÇE - CONTACTO PAPEL

## İÇERİK:

Özel kâğıtlar, balmumu ve nebati yağlar, pigmentler, parafin ve Vazelin.

## KULLANIM TEKNİĞİ:

Oklüzal kuru yüzeye kontakt noktaları kaydetmek için karbon kâğıdı yerleştirin.

## ÖNLEMLER:

Kuru ve iyi havalandırılan bir yerde saklayın.